

EQUIPMENT SPECIFIC TRAINING

(Lateral Cephalostat)

Clinic Name	
Equipment Name	
Operator	
Assessor	

KNOWLEDGE	Tick
Demonstrate how to switch the unit on and perform the warm-up procedure.	
Enter patient details / select from work-list.	
Position patient correctly.	
Use collimation and filters.	
Understand the importance of the use of filters.	
Select correct programme / Select the correct exposure.	
Stamp film in printer with patient's name or write patients name on a label	
Develop film / process the image	
Refill the cassette (if using film)	
Understands the control console and all its functions.	
Find the Protocols, Local Rules and Exposure Charts.	
Demonstrates an ability to use the dose-meter and understands the units of measurement.	
If a failed message comes up when exposing, describe your next action.	
Have knowledge of how to remove the room from action.	
Demonstrates an understanding of infection control.	
Find the emergency stop button & how to reset.	

I certify that the above person has been trained and is competent to use the equipment without supervision. I further certify that this list is not exhaustive and does not constitute comprehensive training.

Signed

Assessor

Date

I have received adequate training and I feel that I am competent to use this equipment without supervision.

Signed

Operator

Date